# ANALITSKE KEMIKALIJE BREZ NEJASNOSTI

### **1 Oznake, standardi in opisi čistosti pri analitskih kemikalijah**

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**Sodium chloride**

puriss. p.a., ACS reagent, reag. ISO, reag. Ph. Eur., ≥99.5%

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**Merck** KGaA Darmstadt, Germany  
Vibrant science and technology

**Merck/MilliporeSigma**

Portfolio brands*	<b>Sigma-Aldrich®</b>	<b>Supelco®</b>	<b>Milli-Q®</b>	<b>Millipore®</b>	<b>SAFC®</b>	<b>BioReliance®</b>
Category	Lab & Production Materials	Analytical Products	Lab Water Solutions	Preparation, Separation, Filtration & Testing Products	Pharma & Biopharma Raw Material Solutions	Pharma & Biopharma Manufacturing Services
Product/ trade name	e.g. Cell-Marque™, Calbiochem®, Isotec®, Msynth®, Trizma®	e.g. Cerilliant®, Certipur®, LiChrosolv®, Spectroquant®, TraceCERT®	e.g. Milli-Q® IQ7000	e.g. Millex®, Mobius®, Amicon®, Clarisolve®, Durapore®, Fractogel®, HY-LiTE®, Hycon®, Snap Id®, Stericup®	e.g. EX-CELL®, CHOZN®, Kilo-Lab®, Parteck®, Sentinel®, Emprove®	e.g. E2E Solutions

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**EMPARTA**

**EMPLURA**

<b>Regulations</b>	ACS, ISO (6353), Reag. Ph Eur	ACS	None
<b>Regulatory support</b>	<ul style="list-style-type: none"> <li>• CoA's</li> <li>• MSDS</li> <li>• BSE/TSE certificates</li> <li>• EMSURE® Quality Documentation</li> </ul>	<ul style="list-style-type: none"> <li>• CoA's</li> <li>• MSDS</li> <li>• BSE/TSE certificates</li> </ul>	<ul style="list-style-type: none"> <li>• CoA's</li> <li>• MSDS</li> </ul>
<b>Purity</b>	99,7 – 99,9%	99,0 – 99,5%	99%
<b>Number of specified parameters</b>	Up to 70	~ 10	4-5
<b>Customer Segments</b>	All regulated industries <ul style="list-style-type: none"> <li>• Pharma</li> </ul>	Less regulated industries <ul style="list-style-type: none"> <li>• Chemical</li> <li>• Food</li> </ul>	Non regulated industries <ul style="list-style-type: none"> <li>• Research</li> <li>• Academia</li> </ul>
<b>Example Applications</b>	<b>Highly demanding lab applications</b> (e.g. Pharma QC)	<b>Routine lab applications</b>	<b>Basic lab applications</b> (e.g. Cleaning, Extraction, synthesis, purification) Chemical Production

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
### Slogan:

*"When you want to be more than sure, choose EMSURE"*

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Document	EMSURE® Premium grade	EMPARTA® Standard grade	EMPLURA® Basic grade
MSDS			
Specification			
Certificate of analysis			
BSE/TSE Certificate			
Enhanced documentation (on request)*			

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# Certificate of Analysis

1.00014.0000 Acetone for analysis EMSURE® ACS,ISO,Reag. Ph Eur  
Batch I1464014

	Spec. Values		Batch Values	
Assay (GC)	≥ 99.8	%	99.9	%
Identity (IR)	conforms		conforms	
Identity B (Ph Eur)	conforms		conforms	
Identity C (Ph Eur)	conforms		conforms	
Appearance of solution	conforms		conforms	
Solubility in water	conforms		conforms	
Water-insoluble substances	conforms		conforms	
Color	≤ 10	Hazen	< 5	Hazen
Density (d 20 °C/20 °C)	0.790 - 0.793		0.791	
Acidity or alkalinity	conforms		conforms	
Titration acid	≤ 0.0002	meq/g	0.0001	meq/g
Titration base	≤ 0.0002	meq/g	< 0.0002	meq/g
Related substances (GC) (Methanol (Impurity A))	≤ 0.05	%	0.03	%
Related substances (GC) (2-Propanol (Impurity B))	≤ 0.05	%	< 0.01	%
Related substances (GC) (Benzene (Impurity C))	≤ 2	ppm	< 1	ppm
Related substances (GC) (other impurity)	≤ 0.05	%	< 0.05	%
Cyclohexane (residual solvent) (GC)	≤ 0.01	%	< 0.01	%
Diacetone (GC)	≤ 0.02	%	< 0.01	%
Ethanol (GC)	≤ 0.01	%	< 0.01	%
Aldehydes (as formaldehyde)	≤ 0.001	%	≤ 0.001	%
Substances reducing potassium permanganate (as O)	≤ 0.00025	%	≤ 0.00025	%
Reducing substances	conforms		conforms	
Chloride (Cl)	≤ 100	ppb	< 20	ppb
Nitrate (NO <sub>3</sub> )	≤ 100	ppb	< 20	ppb
Phosphate (PO <sub>4</sub> )	≤ 100	ppb	< 20	ppb
Sulfate (SO <sub>4</sub> )	≤ 100	ppb	20	ppb
Ag (Silver)	≤ 0.000002	%	≤ 0.000002	%
Al (Aluminium)	≤ 0.00005	%	≤ 0.00005	%
As (Arsenic)	≤ 0.000002	%	≤ 0.000002	%
Au (Gold)	≤ 0.000002	%	≤ 0.000002	%
B (Boron)	≤ 0.000002	%	≤ 0.000002	%
Ba (Barium)	≤ 0.00001	%	≤ 0.00001	%
Be (Beryllium)	≤ 0.000002	%	≤ 0.000002	%
Bi (Bismuth)	≤ 0.000002	%	≤ 0.000002	%
Ca (Calcium)	≤ 0.00005	%	≤ 0.00005	%
Cd (Cadmium)	≤ 0.000005	%	≤ 0.000005	%
Co (Cobalt)	≤ 0.000002	%	≤ 0.000002	%
Cr (Chromium)	≤ 0.000002	%	≤ 0.000002	%

# Certificate of Analysis

1.07021.0000 Acetone for analysis EMPARTA® ACS  
Batch I1464021

	Batch Values	
Purity (GC)	≥ 99.5	%
Identity (IR)	conforms	
Solubility in water	conforms	
Color	≤ 10	Hazen
Titration acid	≤ 0.0003	meq/g
Titration base	≤ 0.0006	meq/g
Methanol (GC)	≤ 0.05	%
2-Propanol (GC)	≤ 0.05	%
Aldehydes (as formaldehyde)	≤ 0.002	%
Substances reducing potassium permanganate (as O)	≤ 0.0003	%
Evaporation residue	≤ 0.001	%
Water	≤ 0.5	%

Date of release (DD.MM.YYYY) 04.12.2025  
Minimum shelf life (DD.MM.YYYY) 30.11.2030

Jeannette David  
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.

# Certificate of Analysis

8.22251.0000 Acetone EMPLURA®  
Batch I1464051

	Spec. Values		Batch Values	
Purity (GC)	≥ 99.0	%	99.9	%
Density (d 20 °C/4 °C)	0.789 - 0.791		0.790	
Evaporation residue	≤ 0.0040	%	< 0.0005	%
Water	≤ 0.30	%	0.03	%

Date of release (DD.MM.YYYY) 04.12.2025  
Minimum shelf life (DD.MM.YYYY) 30.11.2030

Jeannette David  
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.

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Beilstein No.: 3534976

[All Photos](#) (36)

Compare	Product No.	Description	SDS	Pricing
<input type="checkbox"/>	<a href="#">746398</a>	anhydrous, Redi-Dri™, free-flowing, ACS reagent, ≥99%	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">S9888</a>	ACS reagent, ≥99.0%	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">S9625</a>	<i>ReagentPlus</i> ®, ≥99%	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">S3014</a>	Molecular Biology, DNase, RNase, and protease, none detected, ≥99% (titration)	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">S7653</a>	BioXtra, ≥99.5% (AT)	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">S5886</a>	BioReagent, suitable for cell culture, suitable for insect cell culture, suitable for plant cell culture, ≥99%	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">793566</a>	anhydrous, free-flowing, Redi-Dri™, <i>ReagentPlus</i> ®, ≥99%	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">71376</a>	BioUltra, Molecular Biology, ≥99.5% (AT)	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">S1679</a>	meets analytical specification of Ph. Eur., BP, USP, 99.0-100.5%	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">204439</a>	99.999% trace metals basis	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">31434-M</a>	puriss. p.a., ACS reagent, reag. ISO, reag. Ph. Eur., ≥99.5%	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">S6191</a>	BioPerformance Certified, ≥99% (titration), suitable for insect cell culture, suitable for plant cell culture	<a href="#">↓</a>	<a href="#">Expand</a> ▾
<input type="checkbox"/>	<a href="#">71383</a>	tested according to Ph. Eur.	<a href="#">↓</a>	<a href="#">Expand</a> ▾

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## Research Portfolio

### ACS Reagents

**Standard or Basic Grade, Price/Performance**

For **Research** Applications (incl. in process testing)

- High quality for standard lab application
- **Specified acc. to ACS**
- Reduced offer in pack sizes and materials
- In line with nearly all of our competitor's "analytical grade"
- **Specialties like Deuterated and Anhydrous Solvents**

**Supelco**<sup>®</sup>

Analytical Products

## QC Portfolio

### e.g. EMSURE<sup>®</sup>

**Premium Grade**

For **demanding** and **regulated** analytical applications

- Highest quality and reliability
- Specified according to the most important international norms
- Huge variety of pack sizes and packaging materials
- Premium grade for all regulated and demanding lab applications

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


The image displays three overlapping screenshots of the Pregl Digital App interface:

- Dashboard:** Shows a 'Stock' section with four categories: 'Chemicals below minimum' (1), 'Chemicals close to expiration date' (0), 'Chemicals exceeding expiration date' (3), and 'Labware below minimum' (0).
- Drawer Contents:** A list of items with columns for quantity and expiration date. A red header 'Abnormal Conditions A018' is visible.
 

Item	Quantity	Expiration Date
0 Potassium phosphate monoba...	900 g	2021-0
1 Starch	50 g	2024-0
2 Tris	940 g	2022-0
3 SACHAROSE	1000 g	2022-0
4 ETHYLENEDIAMINETETRAA...	428 g	2024-0
5 MAGNESIUM SULFATE HEP...	477 g	2026-0
- Add Alarm:** A configuration screen for setting an alarm. It includes options for 'Alarm type' (Expiration date, Quantity percentage), 'Alarm name' (Half), 'Description' (Half of the stock), and a 'Value in percentage' slider set to 50%.

# PREGL DIGITAL OMOGOČA

**KLIK**

# M-Clarity™ Program


Your Guide to Quality and Portfolio Transparency

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[All Photos \(20\)](#)

Compare	Product No.	Description
<input type="checkbox"/>	<b>1.00014</b>	for analysis EMSURE® ACS,ISO,Reag. Ph Eur
<input type="checkbox"/>	<b>8.22251</b>	for preparative purposes, EMPLURA®

Properties	
grade	ACS reagent
Quality Level	<b>300</b>
Agency	reag. ISO reag. Ph. Eur.
vapor density	2 (vs air)
vapor pressure	184 mmHg ( 20 °C)
description	for analytical purposes for extraction

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Quality Standard ISO 9001

Specifications available

Certificate of Quality or Certificate of Analysis available

Release testing is performed using established protocol

Written SOP for process control

Supplier approval process in line with corporate quality programs

Change notification available as an opt-in for individual products.  
Notifiable events differ between Quality Segments

Release testing is performed using established or published protocol

Site quality self-assessment available

Shelf life/expiration date is identified if applicable

Audits can be requested by customer

Product can be added to a Quality Agreement

Analytical method is verified

Analytical method may be shared upon request

Quality declarations as required by regulation or product application

Process control is verified

Supplier approval by paper audit or questionnaire

Original manufacturer disclosure may be requested with signed confidentiality commitment

Controls for subcontracting are established

Primary packaging component control

Original manufacturer disclosure available with signed confidentiality commitment

Analytical method is validated

Process control is validated

Supplier approval by on-site audit for critical suppliers

Shelf life/expiration date is defined by stability study

Risk based approach to controlled conditions for warehouse & shipping

Original manufacturer disclosure available without confidentiality commitment

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